

TECHNICAL DATA SHEET

V-Print splint comfort

Colour	clear	
Flexural strength	42 MPa*	analogous DIN EN ISO 20795-2
Modulus of elasticity	670 MPa*	analogous DIN EN ISO 20795-2
Water absorption	15 µg/mm ³ *	analogous DIN EN ISO 20795-2
Water solubility	< 2.5 µg/mm ³ *	analogous DIN EN ISO 20795-2
Viscosity	1250 mPa*s	internal testing

*Manufacturing information: SolFlex 3D printer / OtoFlash G171
(Different released printers / post-exposure units may show minor deviations.)

The stated measurement values do not represent a specified value within the scope of the VOCO GmbH quality assurance system.