

TECHNICAL DATA SHEET

V-Print dentbase

Colour	pink	
Flexural strength	90 Mpa*	analogous DIN EN ISO 20795-1
Modulus of elasticity	2.4 GPa*	analogous DIN EN ISO 20795-1
Water absorption	24 $\mu\text{g}/\text{mm}^3$ *	analogous DIN EN ISO 20795-1
Water solubility	< 0.1 $\mu\text{g}/\text{mm}^3$ *	analogous DIN EN ISO 20795-1
Viscosity	1700 mPa*s	internal testing

*Manufacturing information: SolFlex 3D printer / OtoFlash G171
(Different released printers / post-exposure units may show minor deviations.)

The stated measurement values do not represent a specified value within the scope of the VOCO GmbH quality assurance system.