

# TECHNICAL DATA SHEET

## V-Print Try-In

Colour	beige	
Flexural strength	85 MPa*	analogous DIN EN ISO 20795-1
Modulus of elasticity	2.5 GPa*	analogous DIN EN ISO 20795-1
Water absorption	17.5 µg/mm <sup>3</sup> *	analogous DIN EN ISO 20795-1
Water solubility	< 0.1 µg/mm <sup>3</sup> *	analogous DIN EN ISO 20795-1
Viscosity	850 MPa*s	internal testing

\*Manufacturing information: SolFlex 3D printer / OtoFlash G171  
(Different released printers / post-exposure units may show minor deviations.)

The stated measurement values do not represent a specified value within the scope of the VOCO GmbH quality assurance system.